



Wireless Freedom

Focus on scanning your patient without wires getting in the way. Clarius Scanners use Wi-Fi direct to connect to the app on your smart device. Internet access is not required.

Secure Storage

Manage patient images conveniently and securely on the HIPPA-compliant Clarius Cloud, which is included with your scanner. You can also store images directly on your smart device.

Durable and Watertight

With no moving parts or cables and a drop-safe magnesium shell, you don't have to worry about using Clarius in demanding clinical settings. They are sterilizable and submersible to 3 feet.



Manage Your Exams Online for Free

Send your ultrasound exam to Clarius Cloud after you're done scanning. Access your exams stored in the Clarius Cloud from anywhere you have an Internet Connection. Make notes, send exams for review, export PDF reports and more.

Ultrasound Telemedicine Made Simple

With Clarius Live, you can send a link to enable the recipient to view your scan on whatever device they're using.



All-inclusive Pricing

- ✓ No subscription required
- ✓ Unlimited users
- ✓ Unlimited cloud storage

Superior Image Quality

- High-quality images that are comparable to cart-based systems
- Imaging settings are optimized based on scanning application
- Optimal results from specialized scanners

Reliable Portability

- Wireless freedom: no cables to trip over or break, and less risk of infections during procedures
- Stable connectivity from our proprietary multi-channel Wi-Fi architecture

Easy to Use

- Point-and-Shoot Ultrasound™ technology removes need to adjust knobs and switches. Artificial intelligence adjusts parameters in real-time.
- Compact & waterproof device that feels like a standard ultrasound scanner



Model	Frequency
C3 HD	2–6 MHz
L7 HD	4–13 MHz
L15 HD	5–15 MHz
C7 HD	3–10 MHz
PA HD	1–4 MHz
EC7 HD	3–10 MHz

Imaging Modes

B-Mode	Yes
M-Mode	Yes
Power Doppler	Yes
Color Doppler	Yes
Elastography	Upgradeable
Needle Enhance	Upgradeable
Pulsed-Wave Doppler	Upgradeable